Forward and inverse models

Localizing sources using DIPFIT

Robert Oostenveld

Donders Institute for Brain, Cognition and Behaviour Nijmegen, The Netherland





DIPFIT: localizing dipoles

- Motivation
- Ingredients
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models



DIPFIT: localizing dipoles

- <u>Motivation</u>
- Ingredients
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models



Motivation

- Why fit dipoles?
- Why measure EEG?
- Why do ICA?
- Get extra information about brain processes
 - Time course of activity ----> EEG
 - Location of activity ----> fMRI





Difference between EEG and fMRI

- EEG measures post-synaptic potentials
 - related to synchronized neuronal input
- fMRI measures BOLD
 - related to energy consumption
- Different characteristics in the time domain
- Different generators
- Timecourse and location






Why EEG: extra information

- Timecourse
 - ERSP
 - ERP
- Topography
 - Scalp distribution
 - Underlying generators







ICA











Independent component analysis

Mixture of Brain source activity











1





Source modelling



inverse problem




Overview

- Motivation
- Ingredients
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models





Neuronal currents








Symmetry, orientation and activation

radial symmetric

random oriented

asynchronously activated

synchronously activated parallel oriented








Motivation for current dipoles

Neurophysiological motivation





Equivalent current dipoles







17

Motivation for current dipoles

- Neurophysiological motivation
- Physical/mathematical motivation
 - Any current distribution can be written as a multipole expansion
 - First term: monopole (must be zero)
 - Second term: dipole
 - Higher order terms: quadrupole, octupole





Motivation for current dipoles

- Neurophysiological motivation
- Physical/mathematical motivation
 - Any current distribution can be written as a multipole expansion
 - First term: monopole (must be zero)
 - Second term: dipole
 - Higher order terms: quadrupole, octupole
- Convenience
 - dipoles can be used as building block in distributed source models




Overview

- Motivation and background
- Forward modeling
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models
 - Spatial filtering



Volume conductor

- electrical properties of tissue
- geometrical description
- spherical model
- realistic shaped model



 \rightarrow Describes how the currents flow, not where they originate from





Volume conductor

- Advantages spherical model
 - mathematically accurate
 - reasonably accurate
 - computationally fast
 - easy to use
- Disadvantages spherical model
 - inacurate, esp. in some regions
 - difficult alignment with anatomy





Volume conductor

- Advantages realistic model
 - accurate solution for EEG
- Disadvantages realistic model
 - more work
 - individual anatomical MRI required
 - computationally slow(er)
 - numerically instable
 - difficult in interindividual comparison


→The pragmatic solution is to use a standard realistic headmodel for EEG





Realistic volume conductor


- Computational methods for volume conduction problem that allow realistic geometries
 - Boundary Element Method (BEM)
 - Finite Element Method (FEM)
- Geometrical description
 - triangles
 - tetraeders/voxels





Volume conductor: BEM

- Boundary Element Method
 - description of geometry by compartments
 - each compartment is
 - homogenous
 - isotropic
 - important tissues
 - skin
 - skull
 - brain
 - (CSF)
 - triangulated surfaces as boundarie
 - surfaces should be closed





Volume conductor: FEM

- Tesselation of 3D volume in tetraeders
- Large number of elements
- Each tetraeder can have its own conductivity
- FEM is most accurate numerical method
- Computationally expensive
- Accurate conductivities are not (well) known






26

Overview

- Motivation and background
- Forward modeling
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models
 - Spatial filtering



EEG volume conduction





Robert Oostenveld – Forward and inverse models – EEGLAB workshop Aspet 2011



28


EEG volume conduction

- Potential difference between electrodes corresponds
 to current flowing through skin
- Only tiny fraction of current passes through skull
- Therefore the model should describe skull and skin as accurately as possible
- Problems with skull
 - Not visible in anatomical MRI
 - Thickness varies
 - Conductivity is not homogeneous
 - Complex geometry at base of skull





Electric current → magnetic field








30

MEG volume conduction

- Measures sum of fields associated with
 - Primary currents
 - Secondary currents !!!



MEG volume conduction

- Only tiny fraction of current passes through the poorly conductive skull
- Therefore skull and skin can be neglected in the MEG model
- Local conductivity around dipole important
 - geometry
 - conductivity





Differences between EEG and MEG

- more blurred
- deep sources
- electrode noise
- reference electrode
- fixed to head
- skull+skin important for modelling

- no radial or deep sources
- environmental noise
- independent sensors
- head can move
- skull+skin not important





Overview


- Motivation
- Forward modeling
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models





Source modelling

forward problem



Inverse methods

- Single and multiple dipole models
 - Minimize error between model and measured potential/field
- Distributed dipole models
 - Perfect fit of model to the measured potential/field
 - Minimize additional constraint on sources
 - LORETA (smoothness)
 - Minimum Norm (L2)
 - Minimum Current (L1)
- Spatial filtering
 - Scan whole brain with single dipole and compute the filter output at every location
 - MUSIC
 - Beamforming (e.g. LCMV, SAM, DICS)



36

Overview

- Motivation and background
- Forward modeling
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models
 - Spatial filtering





37


Single or multiple dipole models


- Manipulate source parameters to minimize error between measured and model data
 - Location of each source
 - Orientation of each source
 - Strength of each source
- Orientation and strength together correspond to the "dipole moment" and can be estimated linearly
- Position is estimated non-linearly
- Source parameter estimation





Parameter estimation







Parameter estimation: model

measured potential

 $V_i = V(\vec{r_i}) + \text{noise}$

forward model for the data


 $Y_i = Y(r_i; \zeta) + \text{noise}$


select "optimal" model

$$\min_{\zeta} \left\{ \sum_{i=1}^{N} \left(Y(r_i; \zeta) - V(r_i) \right)^2 \right\}$$



Select optimal model









Dipole scanning: grid search

- define grid with allowed dipole locations
- compute optimal dipole moment for each location
- compute value of goal-function
- plot value of goal-function on grid

- number of evaluations:
 - single dipole, 1 cm grid: ~4 000
 - single dipole, $\frac{1}{2}$ cm grid:
 - two dipoles, 1cm grid:

~4 000 ~32 000

~16 000 000



Dipole *fitting*: nonlinear search

- start with an initial guess
- evaluate the local derivative of goal-function
- "walk down hill" to the most optimal solution

number of evaluations: ~100





Overview

- Motivation and background
- Forward modeling
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models
 - Spatial filtering





Distributed source model

- Position of the source is not estimated as such
 - Pre-defined grid (3D volume or on cortical sheet)
- Strength is estimated
 - In principle easy to solve, however...
 - More "unknowns" (parameters) than "knowns" (measurements)
 - Infinite number of solutions can explain the data perfectly
 - Additional constraints required
 - Linear estimation problem





Distributed source model

• Linear estimation

$$\vec{\Psi} = q_1 \vec{\Psi}_1 + q_2 \vec{\Psi}_2 + \dots = \begin{bmatrix} \Psi_{1,1} & \Psi_{2,1} & \cdots \\ \Psi_{1,2} & \Psi_{2,2} & \cdots \\ \vdots & \vdots & \ddots \\ \Psi_{1,N} & \Psi_{2,N} & \cdots \end{bmatrix} \cdot \begin{bmatrix} q_1 \\ q_2 \\ \vdots \end{bmatrix} = \mathbf{L} \cdot \vec{q}$$

$$\vec{q} = \mathbf{L}^{-1} \cdot \vec{\Psi}$$



46

Distributed source model

 $V = L \cdot q + Noise$

$$\min_{q} \{ ||V - L \cdot q ||^2 \} = 0 !!!$$

- Regularized linear estimation: $\min_{q} \{ ||V - L \cdot q ||^{2} + \lambda^{2} \cdot ||D \cdot q ||^{2} \}$
- Constrained linear estimation: $\min_{q} \{q^T \cdot W \cdot q\}$ while $||V - L \cdot q||^2 = 0$



Overview

- Motivation and background
- Forward modeling
 - Source model
 - Volume conductor model
 - Analytical (spherical model)
 - Numerical (realistic model)
 - Comparison EEG and MEG
- Inverse modeling
 - Single and multiple dipole fitting
 - Distributed source models
 - Spatial filtering





Spatial filtering

- position of the source as such is not estimated
- scanning the whole brain
 - single dipole as source
 - estimate activity at each grid location
 - that explains a part of the data
 - that supresses other activity
- various methods
 - <u>multiple signal classification (MUSIC)</u>
 - beamforming
 - LCMV, SAM, DICS, ...
- not a distributed source model, but a distributed representation of the single dipole estimate
- unmixing of data into "signal source" and "noise sources" using assumptions on temporal relation between sources





Summary 1

- Forward modelling
 - Required for the interpretation of scalp topographies
 - Interpretation of scalp topography is "source estimation"
 - Mathematical techniques are available that aid in interpreting scalp topographies -> inverse modeling





5(

Summary 2

- Inverse modeling
 - Model assumption for volume conductor
 - Model assumption for source (I.e. dipole)
 - Additional assumptions on source
 - Single point-like source
 - Multiple point-like sources
 - Distributed source
 - Different mathematical solutions
 - Dipole fitting (linear and nonlinear)
 - Linear estimation (regularized)



