# Independent Component Analysis of Electrophysiological Data



#### **Scott Makeig**

Institute for Neural Computation University of California San Diego

29th EEGLAB Workshop

Bois Perché, Aspet, France

June 17-21, 2019

## Blind EEG Source Separation by Independent Component Analysis



Tony Bell, developer of Infomax ICA

ICA can find distinct EEG source activities -- and their 'simple' scalp maps!

#### Independent Component Analysis of Electroencephalographic Data

Sentt Makrig Naval Health Research Center P.O. Rex 85122 San Diego CA 92188-5122 senttilepl.mag.nbre.navy.ntl


Tayy-Ping Jung
Naval Health Research Center and
Computational Neurobiology Lah
The Salk Institute, P.O. Hoz 85800
San Diego, CA 92186-5800
jung@aalk.4du

Anthony J. Bell Computational Newshiology Lab The Salk Institute, P.O. Ince Si800 San Diego, CA 52185-5800 temptically.edu

Tozenen J. Sejnowski.
Howard Hughen Medical Institute and
Computational Newsokio logy Lah.
The Salk Institute, P.O. Roz 85800
San Diego, OA 92186-3600
terryfeasilt, edn

#### Abstract


Hocause of the distance he tween the skull and havin and their different resistivities, electroeneephologrophie (FFG) data collected from any point on the human scalp includes activity generated within a large hain area. This spatial smearing of PEG data by volume conduction does not involve significant time delays, however, suggesting that the Independent Component Analysis (ICA) a) gorithm of Hell and Sejnowski [1] is suitable for performing hind source separation on PPG data. The KOA algorithm reparator the problem of source identification from that of source localination. First results of applying the ICA algorithm to FEG and escent-related potential (FRP) data collected during a mortained auditory detection tank show: (1) ICA training is insensitive to different random seeds. (2) ICA may be used to segregate obvious artifactual FBG components (line and muscle noise, eye movements) from other sources. (2) ICA is capable of isolating overlapping PPG phenomena, including alpha and theta human and spatially-separable FRP components, to reparate ICA charmels. (4) Nonstationarities in EEG and hehavional state can be tracked using ICA via changes in the amount of residual correlation hetween ICA-filtered output channels.




**CSF** 

## **ICA Assumptions**

- Mixing is **linear** at electrodes
- Propagation delays are negligible
- Component locations are fixed
- Component time courses are independent
- # components <= # scalp channels




### Are EEG effective source signals independent?




#### **Properties of EEG Independent Components**

- Maximally Temporally Independent
- Concurrently Active and Spatially Overlapping
- Dipolar Scalp Maps (Delorme et al., 2012)
- Functionally Distinct
- Between-Subject Similarity / Complexity



## ICA vs. PCA





PCA simply decorrelates the outputs using an **orthogonal mixing matrix**.


PCA makes each successive component account for as much **variance** in the data as possible.

ICA makes each component account for as much **temporally independent information** in the data as possible, with no constraints on the mixing matrix.


## PCA lumps — ICA splits!









#### A P300' visual target response at a vertex channel




#### The vertex channel response sums 238 independent sources.




No more than 30% of any scalp channel variance is produced by any one brain source!




#### ICA finds Non-Brain Independent Component (IC) Processes ...




... separates them from the remainder of the data ...




#### ... and also separates cortical brain IC processes






#### **EEG Source Localization**


LORETA = Low-Resolution Electrical Tomography



Sparse
Compact
Smooth
(SCS)


Distributed IC source location estimation using SCALE-optimal head model.





-----

### ICA is a linear data decomposition method




 $Data = W^{-1} * (W * Data)$ 

## Infomax ICA learning approach

How to make the outputs statistical independent?

Minimize their redundancy or mutual information.

Consider the joint entropy of two components,




(Amari)

#### **Historical Remarks**

- Herault & Jutten ("Space or time adaptive signal processing by neural network models", Neural Nets for Computing Meeting, Snowbird, Utah, 1986): Seminal paper
- Bell & Sejnowski (1995): Information maximization (Infomax)
- Makeig, Bell, Jung, Sejnowski (1996); ICA decomposition of EEG
- Amari et al. (1996): Natural gradient learning
- Cardoso (1996): Joint approximate diagonalization (JADE)
- Hyvarinen (1999): (fastICA)
- Lee/Girolami (1999): Mixture model ICA (Extended Infomax)
- Palmer (2006-9): Adaptive mixture ICA (AMICA)

#### Applications of ICA to biomedical signals

- EEG/ERP analysis (Makeig, Bell, Jung & Sejnowski, NIPS 1996)
- fMRI analysis (McKeown et al., 1998)
- Fetal/mother ECG separation (Cardoso, 1998)
- Electrocorticography (ECoG) (Whitmer, 2010)



## **Important Recent Result (2012)**


Those linear decompositions of multi-channel EEG data that find ICs whose time courses are more temporally independent ...

Also find more ICs whose scalp maps are highly 'dipolar' – i.e., ICs compatible with the spatial projection of a single local cortical (or non-brain, 'artifactual') source process – whose location and distribution can be more or less accurately estimated.

The more independent the component time courses ←→
The larger the number of ~dipolar component scalp maps.

**Hypothesis:** Dipolar ICs = Localized cortical source processes

Pairwise Mutual Information Subspaces




ICA transforms multiple correlated channels into maximally mutually independent component subspaces

#### Distinguishing IC types: ICLabel

### → https://labeling.ucsd.edu

Info: https://sccn.ucsd.edu/wiki/ICLabel



The ICLabel plug-in

#### Multi-model ICA: AMICA



Figure 4: The top panel shows the hypnogram, i.e., sleep stages annotated from the EEG record by a sleep expert, of a sleep session from a single subject. Bottom panels show mean probabilities, within each 30-sec sleep scoring interval, of ICA models learned by an 8-model AMICA decomposition applied to the EEG record. Red-shaded regions highlight changes in model probabilities for relevant models during transitions to and periods of deep sleep (N4). Grayshaded regions highlight probability value changes for relevant models during REM sleep.





Figure 6: Scatter plot of window-mean model probabilities for AMICA model clusters A, B, and C (cf. Fig. 5), each point representing mean model probability within a 30-sec data segment from sleep recordings of 7 healthy subjects and 10 patients. Colors represent expert designated sleep-stage labels for the same data segments. Note the distinct deep sleep (N4) pattern and the relative closeness of wake and REM sleep characteristics.

#### Multi-model ICA: AMICA




# Are effective source *locations* similar across subjects?



S. Makeig (2007)

#### **Visual Working Memory**




Sternberg letter memory task




dipoledensity()

#### **Eyes-closed emotion imagination**




>> dipoledensity()

#### Letter twoback with feedback




#### **Auditory novelty oddball**





Auditory oddball plus novel sounds

#### A. Old/new word memory




#### **B.** Visually cued selective response





#### **Visual Selective Attention Task**




15 subjects


31 channels


Westerfield & Townsend

------




# **Example: frontal midline theta cluster**







## Why don't all subjects contribute to every IC cluster?




## Subject differences?




# Subject differences?





# Subject differences?



#### **Properties of EEG Independent Components**

- Maximally Temporally Independent
- Concurrently Active and Spatially Overlapping
- Dipolar Scalp Maps
- Functionally Distinct
- Between-Subject Similarity / Complexity

# **Important Recent Result**


Those linear decompositions of multi-channel EEG data that find ICs whose time courses are more temporally independent ...

Also find more ICs whose scalp maps are highly 'dipolar' – i.e., ICs compatible with the spatial projection of a single local cortical effective source (or non-brain artifact) – whose cortical location can be accurately estimated given a good forward-problem head model.

**More-independent IC time courses** 




Larger # of dipolar IC scalp maps




# Not the End...

### Multi-model ICA: AMICA



# The 'receptive field' of a bipolar EEG channel!



**Scalp EEG channel** 

Its cortical 'receptive field'

#### Multi-model ICA: AMICA

Reaction speed to simulated driving challenge:




Figure 8: The top panel shows reaction speed changes (inverse of reaction times) in response to lane-departure challenges in one simulated driving session. The three bottom panels show the 5-sec smoothed probabilities of the three ICA models learned by a three-model AMICA decomposition of the whole EEG data session before lane-departure events. Correlation coefficients (r) between each model probability time course and reaction speed are indicated. Black arrows in the lower panel mark brief (alert) periods when model M3 was dominate and reaction speed high.

