



Measure Projection Analysis



Nima Bigdely-Shamlo, Tim Mullen, Ozgur Yigit Balkan Swartz Center for Computational Neuroscience INC, UCSD, 2011

Outline

1. Current EEGLAB Workflow 2. STUDY IC Clustering **3. Issues with IC Clustering** 4. Measure Projection method 5. Practicum

(please start copying the content of workshop USB driver to your computer)

Current EEGLAB Workflow

Single Session Analysis



Study IC Clustering

- Assumes there are functionally equivalent ICs across most subjects.
- Assumes these ICs have similar responses to experimental conditions across all measures (ERP, ERSP, ITC...)
- Creates Non-Overlapping partitions: each IC belongs only to one cluster.

Study IC Clustering



Conceptual Problems with Study IC Clustering

 Components may have similar responses for one measure (e.g. ERSP) but not for the other (e.g. ERP).



Conceptual Problems with Study IC Clustering

Clustering boosts evidence by rejecting ICs that are in the same brain area but show different responses. This makes calculating significance values difficult.

How can we make sure that we are not 'imagining clusters'?

Practical problems with current methods of Study IC Clustering

Large parameter space issue: many different clustering solutions can be produced by changing parameters and measure subsets. Which one should we choose?

EGLAB original	:: Select and compute component measures for later clustering - pop_preclust()
clustering has ~12 parameters	Pre-compute measures on which to cluster components from study 'N400STUDY' Select the cluster to refine during sub-clustering (any existing sub-hierarchy will be overwritten) ParentCluster 1 (151 ICs)
	Pre-compute or Load Dims. Norm. Rel WL. * spectra 10 1 Frequency range [H2] 3 25 * ERPs 10 1 Latency range in ms [10 hi] -2100 1995 * dipoles 3 10 Use channel values Absolute values * Scalp maps 10 1 Use channel values P; [3 25], 'cycles', [3 0 5], 'pe * ERSPs 10 1 Time/freq. parameters P; [3 25], 'cycles', [3 0 5], 'pe * Final dimensions 10 Help Help
	Save STUDY to file ///ata/boinniof/4/emo/5subjecte/N400prexilustaturiv Cancel Help Ok



What are the clusters according to location?



What are the clusters according to circle Size ?



What are the clusters according to both circle location and size? The answer highly depends on how much weight is given to each factor (measure).



Alternatively we could find local neighborhoods (on a grid) with significant (unlikely by chance) similarity in circle Size.



We can define a local-average circle size for each grid location and then cluster these values to form Domains.



- Instead of clustering, we assign to each location in the brain a unique EEG response.
- The response at each location is calculated as the weighted sum of IC responses in its neighborhood.
- Weights are assigned by passing the distance between the location and IC dipole through a Gaussian function.
- The std. of this function represent expected error in dipole localization and inter-subject variability.



- Each EEG measure (ERP, ERSP..) is projected separately.
- Only has one (1) parameter: std. of Gaussian (which has a biological meaning).
- Bootstrap (permutation) statistics can be easily and quickly performed for each point in the brain.
- A regular grid is placed in the brain to investigate every area (with ~8 mm spacing).

- Not all projected values are significant.
- Some are weighted means of ICs with very dissimilar responses.
- Only projected values in neighborhoods with convergent responses are significant.
- Convergence can be expressed as the mean of pairwise similarities in a spatial neighborhood.
- The significance of convergence at each location can be calculated with bootstrapping (permutation).

For a neighborhood with a 'fixed' boundary, for each IC pair we can define a membership function:

 $M(IC_i, IC_j) = M(IC_i)M(IC_j)$

Where M(IC) is one (1) if IC is in the neighborhood and zero (0) otherwise. Convergence can then be defined:

 $convergence = \frac{\sum_{i=1}^{n} \sum_{j=1, j <>i}^{n} \sum_{i=1}^{n} M(i,j)}{\sum_{i=1}^{n} \sum_{j=1, j <>i}^{n} M(i,j)}$

 $\sum^{n} \sum^{n} S(i,j)M(i,j)$

Where M is the neighborhood membership matrix and S is the pairwise similarity matrix. This is basically the mean of pairwise IC similarities around a location in the brain.

Now we can extend this concept of convergence to neighborhoods with 'soft' Gaussian boundaries, for each IC pair we modify the membership function:

 $MG(IC_i, IC2_j) = MG(IC_i)MG(IC_j)$ Where $MG(IC) = e^{-\frac{d^2}{2\sigma^2}}$ distance from IC equiv. dipole to neighborhood center). Convergence can now be defined as:



Where S is the pair-wise similarity matrix. This is basically the weighted mean of IC similarities around a location in the brain.



To better visualize measure responses in areas with significant convergence, they can be summarized into different *domains*. The exact number of these domains depends on how similar their exemplars are allowed to be.

Below you can see ERSP responses in an EEG experiment form three (3) domains.





Rapid Serial Visual Presentation Experiment

•8 subjects

15 Sessions

Visual target detection
257 components with equiv. dipoles inside the brain





Clusters Cluster 11, Target ERSP Cluster 14, Target ERSP Cluster 4, Target ERSP 400 Time (ms) 400 Time (ms) -400 800 800 400 -4000 -4000 0 Time (ms) Cluster 1, Target ERSP Cluster 15, Target ERSP Cluster 12, Target ERSP (Hz)

> 400 Time (ms)

800

-400

400 Time (ms) 800

-400

requency (Hz)

13

-400

0 400 Time (ms) 800





Subject Space

Measure or dipole density similarity between each two EEG subjects (or sessions) may be averaged over a region of interest (ROI) and visualize using multi-dimensional scaling.



Projected ERSP at ROIs







Measure Projection: Summary

- Enables us to compare subjects, groups and conditions at every brain location.
- Enables us to calculate significance on every step.
- Enables us to perform new types of analysis that we could not do with IC clusters (e.g. subject similarity space)
- All types of analysis that can be done on IC clusters, can also be performed in Measure Projection framework.